

## Return Packing Instructions v1.00


When returning items to AMTI for repair or recalibration please pack them well so that they are padded sufficiently and do not shift inside the container you are returning them in.

It is best to pack them in the original container they came in along with the packing material used (foam peanuts or foam pieces).

**When returning a force platform please remove the mounting plate or mounting rails before placing the platform in the shipping container.**

If a shipping box is provided by AMTI:

1. Place 4 foam pieces in each corner of the shipping box that AMTI provided.
2. Place the force platform in the box resting it on the foam corner pieces.
3. Place the additional 4 foam pieces on the top of the force platform.
4. Cover the box and tape as shown below:



Note: Each of the 4 locations of the filament tape (shown in blue) should wrap entirely around the box at least two times.